

Background

The association between threshold clozapine blood levels and positive clinical response has been demonstrated in many studies.¹⁻⁷ Practice guidelines, clinical support tools, and expert consensus recommendations⁸⁻¹² recommend measuring clozapine levels for both therapeutic drug monitoring and adherence monitoring. Measurement of clozapine blood levels can provide clinicians objective evidence to help avoid treatment failures. Knowledge of nonadherence to atypical antipsychotics, including clozapine, may help clinical decision making.¹³

Point of Care (POC) testing of clozapine blood levels offers:

- ❖ immediacy of results for clinicians and patients,⁹
- ❖ capillary sampling that is preferred by both clinicians and patients,¹⁴ and
- ❖ improved utilization of clozapine, which is the only effective drug for treatment-resistant schizophrenia, where monitoring levels is important for effective use.¹⁵

The aim of this study was to validate the POC MyCare Insite Clozapine Test for use in clinical practice.

Results

Assay bias associated with endogenous interferents and an extreme range of hematocrits was < 10% (Table 1).

The specificity of the antibody has previously been characterized: no interference from over 140 small molecules was found.¹⁶ The assay is specific for clozapine and does not cross-react with metabolites or structurally related antipsychotic medications (Table 2).

Results are reported from 0 to 1600 ng/mL. Sensitivity (Table 3) was determined according to CLSI Guideline EP17-A2.

The agreement between Insite and LC-MS/MS observed for the calibrator assignment training set (Chart 1) was reproduced in the validation method comparison (Chart 2). Four samples > 1600 ng/mL on the MyCare Insite were eliminated from the analysis, as were five discrepant samples to be retested. CLSI EP6 linearity was demonstrated with 10 levels from 122 to 1390 ng/mL, measured n = 5 (Chart 3). Precision was evaluated over 20 days according to CLSI guideline EP05-A3; the average results from 3 lots generated by 6 different operators using 22 analysers are shown (Chart 4). Repeatability with patient samples (Table 4) showed that replicate two differed from replicate one by no more than 10%.

Table 1: Interferents and Hematocrit

Interferent	Level	
Rheumatoid Factor	510 IU/mL	
Total Protein	6.1 g/dL	61 g/L
Icteric	44 mg/dL	753 µmol/L
Lipemia	510 mg/dL	5.7 mmol/L
Hemolysate	1050 mg/dL	
Hematocrit	35% - 50%	

Table 2: Selectivity

Compound	Cross-Reactivity
Norclozapine	5.5%
Clozapine N-oxide	5.3%
Quetiapine	< 0.5%
Olanzapine	6.6%

Table 3: Sensitivity

Limit of	ng/mL
Blank (LoB)	34
Detection (LoD)	126
Quantitation (LoQ)	175

Chart 1: Method Comparison for Calibrator Assignment, n = 95, R = 0.9, slope = 1.02, intercept = -23




Chart 2: Method Comparison for Validation, n = 304, R = 0.9, slope = 0.971, intercept = -21.2




Chart 3: Linearity

0% deviation from linearity; ≤ 10% CV above LoQ





Table 4: Patient Sample Repeatability

ng/mL (mean)	Difference
453	-1%
382	-10%
329	-2%
199	1%

Chart 4: Precision

repeatability < 12% on whole blood samples; within-laboratory < 16%



Conclusions

The MyCare Insite Clozapine Test:

- Is the first point of care method developed to test patients on clozapine
- Provides accurate and reliable point-of-care results for clozapine
- Could be used to meet an unmet clinical need by providing psychiatrists with close to patient testing and immediate results for clozapine levels.

