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







## Bibliografija

- Janssen Pharmaceuticals, Inc. Risperdal® (risperidone) Prescribing information. 2017.
- Janssen Pharmaceuticals, Inc. Risperdal Consta® (risperidone) Long-acting Injection Prescribing Information. 2017.
- Janssen Pharmaceuticals, Inc. Invega® (paliperidone) Prescribing information. 2017.
- Janssen Pharmaceuticals, Inc. Invega Sustenna® (paliperidone palmitate) extended-release injectable suspension Prescribing Information. 2017.
- Seto K, Dumontet J, Ensom MH. Risperidone in schizophrenia: is there a role for therapeutic drug monitoring? *Ther Drug Monit.* 2011;33(3):275-283.
- Velligan DI, Weiden PJ, Sajatovic M, et al. Assessment of adherence problems in patients with serious and persistent mental illness: recommendations from the Expert Consensus Guidelines. *J Psychiatr Pract.* 2010;16(1):34-45.
- Higashi K, Medic G, Littlewood KJ, Diez T, Granstrom O, De Hert M. Medication adherence in schizophrenia: factors influencing adherence and consequences of nonadherence, a systematic literature review. *Ther Adv Psychopharmacol.* 2013;3(4):200-218.
- Haddad PM, Brain C, Scott J. Nonadherence with antipsychotic medication in schizophrenia: challenges and management strategies. *Patient Relat Outcome Meas.* 2014;5:43-62.
- Hiemke C, Bergemann N, Clement HW, et al. Consensus Guidelines for Therapeutic Drug Monitoring in Neuropsychopharmacology: Update 2017. *Pharmacopsychiatry.* 2018;51:9-62.
- Bowskill SVJ, Handley SA, Fisher DS, Flanagan RJ, Patel MX. Risperidone and Total 9-Hydroxyrisperidone in Relation to Prescribed Dose and Other Factors: Data From a Therapeutic Drug Monitoring Service, 2002–2010. *Ther Drug Monit.* 2012;34(3):349-355.
- Lostia AM, Mazzarini L, Pacchiarotti I, et al. Serum levels of risperidone and its metabolite, 9-hydroxyrisperidone: correlation between drug concentration and clinical response. *Ther Drug Monit.* 2009;31(4):475-481.
- Aravagiri M, Marder SR, Nuechterlein KH, Gitlin MJ. Intra- and interindividual variations in steady-state plasma concentrations of risperidone and 9-hydroxyrisperidone in schizophrenic patients treated chronically with various doses of risperidone. *Ther Drug Monit.* 2003;25(6):657-664.
- Gefvert O, Eriksson B, Persson P, et al. Pharmacokinetics and D2 receptor occupancy of long-acting injectable risperidone (Risperdal Consta) in patients with schizophrenia. *Int J Neuropsychopharmacol.* 2005;8(1):27-36.
- Melkote R, Singh A, Vermeulen A, Remmerie BM, Savitz A. Relationship between Antipsychotic Blood Levels and Treatment Failure during the CATIE Study. International Congress on Schizophrenia Research. 2017;Poster.
- CLSI. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute, 2014.
- CLSI. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- Mannens G, Huang ML, Meuldermans W, Hendrickx J, Woestenborghs R, Heykants J. Absorption, metabolism, and excretion of risperidone in humans. *Drug Metab Dispos.* 1993;21(6):1134-1141.
- Spina E, Cavallaro R. The pharmacology and safety of paliperidone extended release in the treatment of schizophrenia. *Expert Opin Drug Saf.* 2007;6(6):651-662.
- Vermeir M, Naessens I, Remmerie B, et al. Absorption, metabolism, and excretion of paliperidone, a new monoaminergic antagonist, in humans. *Drug Metab Dispos.* 2008;36(4):769-779.
- CLSI. Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition CLSI document EP7-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
- NCCLS. Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. NCCLS document EP6-A. Wayne, PA: NCCLS; 2003.
- CLSI. Measurement Procedure and Bias Estimation Using Patient Samples; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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EMERGO EUROPE  
Westervoortsedijk 60,  
6827 AT Arnhem  
The Netherlands

Casus Switzerland GmbH  
Hinterbergstrasse 49  
6312 Steinhausen  
Switzerland

Saladax Biomedical, Inc.  
116 Research Dr.  
Bethlehem, PA 18015 SAD  
[www.saladax.com/bci\\_applications/](http://www.saladax.com/bci_applications/)

Odgovorna osoba u UK-u:  
Emergo Consulting (UK) Limited  
c/o Cr360 – UL International  
Compass House, Vision Park Histon  
Cambridge CB24 9BZ  
Ujedinjena Kraljevina

Naručitelj iz Australije  
ACRA Regulatory Services Pty Ltd  
7/ 84 Poinciana Avenue,  
Tewantin, QLD 4565 Australija

Naručitelj iz Novog Zelanda  
ACRA Regulatory Services Limited  
182 Teasdale Street,  
Te Awamutu, 3800, Novi Zeland

© 2023, Saladax Biomedical, Inc.